

Plant Culture Bottle II (#310501)

1. Description

Plant Culture Bottles are used for culture of sapling or callus.

2. Features

- ❖ Durable plastic bottle
- ❖ Designed for reduced contamination
- ❖ Screw cap
- ❖ Autoclavable

3. Product Range

Cat. No.	Material (Bottle / Cap)	Dimension (mm)	Total Volume (ml)	Sterile	Packaging
310501	PC / HDPE	115.00 x115.00 x 140.00	500.00	-	1 / 50

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs



SPL-TDS-PCB2


For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 21th, 2016

6. Specifications & Technical Drawing

- Specifications

Dimensions		Material (Resin)	Packaging
Lid	Dish		
Diameter Inner:	57.00 Ø	Crystal-grade polycarbonate (PC)	<ul style="list-style-type: none"> ▪ Pieces/Sleeve: 1 ▪ Pieces/Case: 50
Height:	17.00 mm		
	140.00 mm		

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-PCB2

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 21st, 2016