

TECHNICAL DATA SHEET

Part Number	33350-100G-HON
Description	EDTA dipotassium magnesium salt

Technical specifications :

- ✓ ethylenediaminetetraacetic acid dipotassium magnesium salt =98% 100g
- ✓ CAS N°15708-48-2
- ✓ Chemical purity ≥98%
- ✓ Molecular weight 426.74 g/mol
- ✓ Molecular formula $C_{10}H_{12}K_2MgN_2O_8 \cdot 2H_2O$
- ✓ Flash point Not applicable
- ✓ Material form solid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 100g

- ✓ P Statements P280

Sales unit :

100g EDTA dipotassium magnesium salt
