

TECHNICAL DATA SHEET

Part Number	33545-1L-HON
Description	Buffer solution pH 6.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 6,0, citric acid / NaOH sol., with fungicide, NIST traceable 1l
- ✓ CAS N° 100-56-1
- ✓ Density 1.00 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1L

- ✓ P Statements P280

Sales unit :

1L Buffer solution pH 6.0 (20 °C)
