

TECHNICAL DATA SHEET

Part Number	33552-1L-HON
Description	Buffer solution pH 13.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 13,0, glycine / sodium hydroxide / sodium chloride sol. 1l
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS05
- ✓ H Statements H290
- ✓ P Statements P280 ; P234

Sales unit :

1L Buffer solution pH 13.0 (20 °C)
