

TECHNICAL DATA SHEET

Part Number	33582-1L-HON
Description	Ammonia buffer solution

Technical specifications :

- ✓ ammonia buffer sol., ammonium chloride / ammonia, ph 10, for complexometry 1l
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS07 ; GHS05
- ✓ H Statements H412 ; H314 ; H335
- ✓ P Statements P284 ; P280 ; P260

Sales unit :

1L Ammonia buffer solution
