

TECHNICAL DATA SHEET

Part Number	33595-1L-HON
Description	Buffer solution pH 9.180 (25°C)

Technical specifications :

- ✓ buffer solution pH 9,180 (25°C), with fungicide, according to din 19266 1I
- ✓ CAS N°1330-43-4
- ✓ Density 1.0 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1L

- ✓ P Statements P280

Sales unit :

1L Buffer solution pH 9.180 (25°C)
