

TECHNICAL DATA SHEET

| | |
|--------------------|------------------------------------|
| Part Number | 33597-1L-HON |
| Description | Buffer solution pH 7.413 (25°C) |

Technical specifications :

- ✓ buffer solution ph 7,413 (25°c), with fungicide, according to din 19266 1l
- ✓ Density 1.000 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1L
- ✓ P Statements P280

Sales unit :

1L Buffer solution pH 7.413 (25°C)
