

TECHNICAL DATA SHEET

Part Number	33643-100ML-HON
Description	Buffer solution pH 4.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 4,0, with fungicide, citric acid/naoh/nacl sol. nist trac. 100ml
- ✓ Density 1.0 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 100mL

- ✓ P Statements P280

Sales unit :

100mL Buffer solution pH 4.0 (20 °C)
