

TECHNICAL DATA SHEET

Part Number	33643-10L-VP-HON
Description	Buffer solution pH 4.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 4,0, with fungicide, citric acid/naoh/nacl sol., nist trac. 10l
- ✓ Density 1.0 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type VOLPAC™ Container
- ✓ Package size 10L

- ✓ P Statements P280

Sales unit :

10L Buffer solution pH 4.0 (20 °C)
