

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	33643-1L-HON
<b>Description</b>	Buffer solution pH 4.0 (20 °C)

---

## Technical specifications :

- ✓ buffer sol. pH 4,0, with fungicide, citric acid/naoh/nacl sol., nist trac. 1l
- ✓ CAS N°18996-35-5
- ✓ Density 1.0 g/cm<sup>3</sup>
- ✓ Flash point Not applicable
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1L
  
- ✓ P Statements P280

## Sales unit :

1L Buffer solution pH 4.0 (20 °C)

---