

# TECHNICAL DATA SHEET

---

|                    |                                |
|--------------------|--------------------------------|
| <b>Part Number</b> | 33643-5L-VP-HON                |
| <b>Description</b> | Buffer solution pH 4.0 (20 °C) |

---

## Technical specifications :

- ✓ buffer sol. pH 4,0, with fungicide, citric acid/naoh/nacl sol., nist trac. 5l
- ✓ Density 1.0 g/cm<sup>3</sup>
- ✓ Flash point Not applicable
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type VOLPAC™ Container
- ✓ Package size 5L
  
- ✓ P Statements P280

## Sales unit :

5L Buffer solution pH 4.0 (20 °C)

---