

TECHNICAL DATA SHEET

Part Number	33646-10L-VP-HON
Description	Buffer solution pH 7.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 7,0, k₂hpo₄ / na₂hpo₄, nist traceable, with fungicide 10l
- ✓ CAS N°7558-79-4
- ✓ Density 1 g/cm³
- ✓ Flash point does not flash
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type VOLPAC™ Container
- ✓ Package size 10L

- ✓ P Statements P280

Sales unit :

10L Buffer solution pH 7.0 (20 °C)
