

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	33646-500ML-HON
<b>Description</b>	Buffer solution pH 7.0 (20 °C)

---

## Technical specifications :

- ✓ buffer sol. pH 7,0, k<sub>2</sub>hpo<sub>4</sub> / na<sub>2</sub>hpo<sub>4</sub>, nist traceable, with fungicide 500ml
- ✓ CAS N°7558-79-4
- ✓ Density 1 g/cm<sup>3</sup>
- ✓ Flash point does not flash
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 500mL
  
- ✓ P Statements P280

## Sales unit :

500mL Buffer solution pH 7.0 (20 °C)

---