

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	33648-500ML-HON
<b>Description</b>	Buffer solution pH 9.0 (20 °C)

---

## Technical specifications :

- ✓ buffer solution pH 9,0, borax / hydrochloric acid, nist traceable 500ml
- ✓ Density 1 g/cm<sup>3</sup>
- ✓ Flash point Not applicable
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 500mL
  
- ✓ P Statements P280

## Sales unit :

500mL Buffer solution pH 9.0 (20 °C)

---