

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	33649-500ML-HON
<b>Description</b>	Buffer solution pH 10.0 (20°C)

---

## Technical specifications :

- ✓ buffer solution ph 10,0 (20°c), borax / sodium hydroxide 500ml
- ✓ Density 1.01 g/cm<sup>3</sup>
- ✓ Flash point Not applicable
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 500mL
  
- ✓ P Statements P280

## Sales unit :

500mL Buffer solution pH 10.0 (20°C)

---