

TECHNICAL DATA SHEET

Part Number	33650-1L-HON
Description	Buffer solution pH 11.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 11,0, boric acid / sodium hydroxide sol. / potassium chloride 1l
- ✓ CAS N°13840-56-7
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid, clear

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1L

- ✓ P Statements P280

Sales unit :

1L Buffer solution pH 11.0 (20 °C)
