

TECHNICAL DATA SHEET

Part Number	33650-500ML-HON
Description	Buffer solution pH 11.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 11.0, h3bo3 / naoh sol. / kcl 500ml
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid, clear

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 500mL

- ✓ P Statements P280

Sales unit :

500mL Buffer solution pH 11.0 (20 °C)
