

TECHNICAL DATA SHEET

| | |
|--------------------|---------------------------------|
| Part Number | 33651-100ML-HON |
| Description | Buffer solution pH 12.0 (20 °C) |

Technical specifications :

- ✓ buffer sol. pH 12.0, di-sodium hydrogen phosphate / sodium hydroxide sol. 100ml
- ✓ CAS N°7601-54-9
- ✓ Density 1.0 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 100mL

- ✓ P Statements P280

Sales unit :

100mL Buffer solution pH 12.0 (20 °C)
