

TECHNICAL DATA SHEET

Part Number	33665-1L-HON
Description	Buffer solution pH 4.0 (20 °C)

Technical specifications :

- ✓ buffer solution ph 4,0 red colored with fungicide 1l
- ✓ CAS N°100-56-1
- ✓ Density 1.0 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 6 x 1L

- ✓ P Statements P280

Sales unit :

6 x 1L Buffer solution pH 4.0 (20 °C)
