

TECHNICAL DATA SHEET

Part Number	33665-500ML-HON
Description	Buffer solution pH 4.0 (20 °C)

Technical specifications :

- ✓ buffer solution ph 4,0 red colored with fungicide 500ml
- ✓ Density 1.0 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 6 x 500mL

- ✓ P Statements P280

Sales unit :

6 x 500mL Buffer solution pH 4.0 (20 °C)
