

TECHNICAL DATA SHEET

Part Number	33665-5L-VP-HON
Description	Buffer solution pH 4.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 4,0, red colored, citric acid/naoh/nacl sol., with fungicide 5l
- ✓ Density 1.0 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type VOLPAC™ Container
- ✓ Package size 5L

- ✓ P Statements P280

Sales unit :

5L Buffer solution pH 4.0 (20 °C)
