

TECHNICAL DATA SHEET

Part Number	33666-5L-VP-HON
Description	Buffer solution pH 7.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 7,0, green colored, K_2HPO_4 / Na_2HPO_4 , with fungicide 5l
- ✓ Density 1.0 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type VOLPAC™ Container
- ✓ Package size 5L

- ✓ P Statements P280

Sales unit :

5L Buffer solution pH 7.0 (20 °C)
