

TECHNICAL DATA SHEET

Part Number	33668-1L-HON
Description	Buffer solution pH 10.0 (20°C)

Technical specifications :

- ✓ buffer solution pH 10,0 (20°C), violet colored, borax / sodium hydroxide 1l
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 6 x 1L

- ✓ P Statements P280

Sales unit :

6 x 1L Buffer solution pH 10.0 (20°C)
