

TECHNICAL DATA SHEET

Part Number	34607-500ML-HON
Description	Phenolphthalein solution

Technical specifications :

- ✓ phenolphthalein solution, indicator, reag. ph. eur., 1% in ethanol 500ml
- ✓ Molecular formula C₂₀H₁₄O₄
- ✓ no data available
- ✓ Flash point 14 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 500mL

- ✓ GHS pictograms GHS07 ; GHS08 ; GHS02
- ✓ H Statements H225 ; H319 ; H341 ; H350
- ✓ P Statements P280 ; P201 ; P210 ; P243

Sales unit :

500mL Phenolphthalein solution
