

Standard Analysis Certificate

PRODUCT :GLYCEROL 90% (28°Bé) RE PURO
CODE :346130
METHOD :8806

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Appearance of solution	-	Conform
Acidity or alkalinity	-	Conform
Esters	-	Conform
Sugar	-	Conform
Density at 20° C	-	1.220 ÷ 1.233
Refractive index at 20°C	-	1.4490 ÷ 1.4550
Water (K.F.)	%	11.5 ÷ 16.5
Aldehydes	ppm	<= 10
Chloride	ppm	<= 10
Chlorinated compounds	ppm	<= 30
Heavy metals (Pb)	ppm	<= 5
Residue on ignition	ppm	<= 100
Assay	%	83.5 ÷ 88.5

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)