

TECHNICAL DATA SHEET

Part Number	34629-2.5L-HON
Description	Silver sulfate-Sulfuric acid solution

Technical specifications :

- ✓ silver sulfate sulfuric acid sol., volum. sol. for cod det. din 38409 2,5l
- ✓ CAS N°7664-93-9
- ✓ Density 1.85 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 2.5L

- ✓ GHS pictograms GHS05 ; GHS09
- ✓ H Statements H410 ; H314
- ✓ P Statements P260 ; P280 ; P284

Sales unit :

2.5L Silver sulfate-Sulfuric acid solution
