

Honeywell

Honeywell Specialty Chemicals Seelze GmbH
Wunstorferstrasse 40
Seelze, 30926
Telefon: +49 5137 999-0
lab.honeywell.com

Brand: Honeywell Fluka
Product: 34629

CAS Numbers (All)
7664-93-9

Silver sulfate solution 10 g/l in sulfuric acid for COD determination DIN38409, part 41

Parameter	Specification	Units
assay of H ₂ SO ₄	94-96	%
assay of Ag ₂ SO ₄	0.54-0.56	%
density (20°C)	1.83-1.85	g/ml

This document was produced electronically and is valid without a signature