

TECHNICAL DATA SHEET

Part Number	34698-1L-HON
Description	HYDRANAL™ - Medium K

Technical specifications :

- ✓ hydranal medium k, medium 1-comp. k. fischer titr. in -cho and ket. w/o meoh 1l
- ✓ CAS N°101023-55-6
- ✓ Density 1.255 g/cm³ (20 °C)
- ✓ Flash point 56 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS05 ; GHS08 ; GHS06 ; GHS02
- ✓ H Statements H372 ; H360Fd ; H351 ; H226 ; H302 ; H315 ; H318 ; H331 ; H336
- ✓ P Statements P280 ; P260 ; P284

Sales unit :

1L HYDRANAL™ - Medium K
