

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	34726-500ML-HON
<b>Description</b>	HYDRANAL™ - Coulomat E

---

## Technical specifications :

- ✓ hydranal coulomat e analyte for kf (etoh) for cells w/ or w/o diaphr. 500ml
- ✓ Density 0.980 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point 12 °C
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 500mL
  
- ✓ GHS pictograms GHS07 ; GHS05 ; GHS08 ; GHS02
- ✓ H Statements H225 ; H302 ; H332 ; H318 ; H360D ; H371 ; H373
- ✓ P Statements P284 ; P260 ; P280

## Sales unit :

500mL HYDRANAL™ - Coulomat E

---