

TECHNICAL DATA SHEET

Part Number	34734-1L-HON
Description	HYDRANAL™ - CompoSolver E

Technical specifications :

- ✓ hydranal composolver e, medium for volum. 1-comp. k. fischer titr. (etoh) 1l
- ✓ Density 0.802 g/cm³ (20 °C)
- ✓ Flash point 12 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS08 ; GHS07 ; GHS02
- ✓ H Statements H319 ; H225 ; H351 ; H360D
- ✓ P Statements P284 ; P280

Sales unit :

1L HYDRANAL™ - CompoSolver E
