

TECHNICAL DATA SHEET

Part Number	34741-1L-HON
Description	HYDRANAL™ - Methanol dry

Technical specifications :

- ✓ hydranal methanol dry, medium for volumetric 1-comp. k. fischer titr. 1l
- ✓ CAS N°67-56-1
- ✓ Molecular weight 32.04 g/mol
- ✓ Density 0.79 g/cm³ (20 °C)
- ✓ Flash point 11 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS02 ; GHS06 ; GHS08
- ✓ H Statements H370 ; H301 ; H311 ; H331 ; H225
- ✓ P Statements P260 ; P280 ; P284

Sales unit :

1L HYDRANAL™ - Methanol dry
