

TECHNICAL DATA SHEET

Part Number	34800-1L-HON
Description	HYDRANAL™ - Solvent

Technical specifications :

- ✓ hydranal solvent, medium for volumetric 2-comp. k. fischer titr. (meoh based) 1l
- ✓ Density 0.890 g/cm³ (20 °C)
- ✓ Flash point 11 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS05 ; GHS08 ; GHS06 ; GHS02
- ✓ H Statements H370 ; H360D ; H225 ; H301 ; H311 ; H331 ; H314
- ✓ P Statements P284 ; P280 ; P260

Sales unit :

1L HYDRANAL™ - Solvent
