

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	34801-1L-HON
<b>Description</b>	HYDRANAL™ - Titrant 5

---

## Technical specifications :

- ✓ hydranal titrant 5for volum. 2-comp. k. fischer titr. (meoh based) 1l
- ✓ CAS N°7553-56-2
- ✓ Density 0.850 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point 11 °C
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 1L
  
- ✓ GHS pictograms GHS08 ; GHS06 ; GHS02
- ✓ H Statements H373 ; H370 ; H301 ; H311 ; H331 ; H225
- ✓ P Statements P260 ; P280 ; P284

## Sales unit :

1L HYDRANAL™ - Titrant 5

---