

TECHNICAL DATA SHEET

Part Number	34805-1L-HON
Description	HYDRANAL™ - Composite 5

Technical specifications :

- ✓ hydranal composite 5, titr. agent 1-comp. k. fischer titr. (meoh free) 1l
- ✓ Density 1.17 g/cm³ (25 °C)
- ✓ Flash point 99 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS08
- ✓ H Statements H373 ; H360D ; H351
- ✓ P Statements P260 ; P280 ; P284

Sales unit :

1L HYDRANAL™ - Composite 5
