

TECHNICAL DATA SHEET

Part Number	34806-1L-HON
Description	HYDRANAL™ - Composite 2

Technical specifications :

- ✓ hydranal composite 2, titr. agent 1-comp. k. fischer titr. (meoh free) 1l
- ✓ CAS N°68007-08-9
- ✓ Density 1.110 g/cm³ (20 °C)
- ✓ Flash point 103 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS08
- ✓ H Statements H351 ; H360D ; H373
- ✓ P Statements P260 ; P280 ; P284

Sales unit :

1L HYDRANAL™ - Composite 2
