

# TECHNICAL DATA SHEET

---

|                    |                        |
|--------------------|------------------------|
| <b>Part Number</b> | 34807-500ML-HON        |
| <b>Description</b> | HYDRANAL™ - Coulomat A |

---

## Technical specifications :

- ✓ hydranal coulomat a anolyte for kf (methanol-chloroform-based) 500ml
- ✓ CAS N°68007-08-9
- ✓ Density 1.070 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point 11 °C
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 500mL
  
- ✓ GHS pictograms GHS02 ; GHS06 ; GHS08 ; GHS05
- ✓ H Statements H372 ; H370 ; H360D ; H225 ; H301 ; H311 ; H331 ; H314 ; H351
- ✓ P Statements P284 ; P280 ; P260

## Sales unit :

500mL HYDRANAL™ - Coulomat A

---