

TECHNICAL DATA SHEET

Part Number	34811-1L-HON
Description	HYDRANAL™ - Titrant 2

Technical specifications :

- ✓ hydranal titrant 2 for volum. 2-comp. k. fischer titr. (meoh based) 1l
- ✓ Density 0.800 g/cm³ (20 °C)
- ✓ Flash point 11 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS08 ; GHS06 ; GHS02
- ✓ H Statements H373 ; H370 ; H225 ; H301 ; H311 ; H331
- ✓ P Statements P260 ; P280 ; P284

Sales unit :

1L HYDRANAL™ - Titrant 2
