

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	34816-1L-HON
<b>Description</b>	HYDRANAL™ - Composite 5 K

---

## Technical specifications :

- ✓ hydranal composite 5 k titrating agent 1-component titr. aldehydes ketones 1l
- ✓ Density 1.170 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point 99 °C
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 1L
  
- ✓ GHS pictograms GHS08
- ✓ H Statements H351 ; H360D ; H373
- ✓ P Statements P284 ; P280 ; P260

## Sales unit :

1L HYDRANAL™ - Composite 5 K

---