

TECHNICAL DATA SHEET

Part Number	34817-1L-HON
Description	HYDRANAL™ - Working Medium K

Technical specifications :

- ✓ hydranal working medium k one-component titration in aldehydes and ketones
1l
- ✓ CAS N°107-07-3
- ✓ Density 1.350 g/cm³ (20 °C)
- ✓ Flash point 100 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS06 ; GHS08 ; GHS05 ; GHS09
- ✓ H Statements H411 ; H372 ; H361d ; H351 ; H336 ; H318 ; H315 ; H300 ; H310 ; H330
- ✓ P Statements P284 ; P260 ; P280

Sales unit :

1L HYDRANAL™ - Working Medium K
