



PRODUCT CODE: 252531

Gram-Hucker's Safranin O solution for clinical diagnosis

Sanitary product for in vitro diagnostic class A.

TARIC 3204 13 00 90

PHYSICAL DATA: liquid, Clear, cherry red, Miscible with alcohol and water • D 20/4 0,99 •

BIBLIOGRAPHY: Clark G., "Staining Procedures", 4th Ed., Williams and Wilkins ,Baltimore , 377 (1981) ; R. García del Moral , "Laboratorio de AnatomíaPatológica", 1ª Ed., McGraw-Hill , 136 (1993) •

HAZARDOUS: VLA-ED (ethanol) 1.000 ppm1.910 mg/m3

WEIGHT/VOLUME INFORMATION: 1l~0,990 kg 1kg~1,01 l

COMPOSITION:

Safranin O	0,25 g
Ethanol absolute	10 ml
Water (s.q.m)	100 ml

SPECIFICATIONS:

Maximum limit of impurities

ABS	
at λ at 519 \pm 3 nm; dil. 1:1000	0,200 - 0,285
Staining suitability	passes test

