



PRODUCT CODE: 252164

Buffer Solution pH 7.2 for clinical diagnosis
for hematology, according to Weise

TARIC 3822 00 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,001 • pH7,2 •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,001 kg 1kg~0,999 l

COMPOSITION:

Potassium di-Hydrogen Phosphate	40 mg
di-Sodium Hydrogen Phosphate 12-hydrate	151 mg
Water (s.q.m)	100 ml

SPECIFICATIONS:

Identity :	
Identity	passes test
pH	7,0 - 7,4