

TECHNICAL DATA SHEET

Part Number	34828-40ML-HON
Description	HYDRANAL™ - Water Standard 1.0

Technical specifications :

- ✓ hydranal water standard 1.0 for titration water content 1 mg/g = 0,1%, 4ml x10
- ✓ Density 1.000 g/cm³ (20 °C)
- ✓ Flash point 43 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Ampoule
- ✓ Package size 40mL

- ✓ GHS pictograms GHS02
- ✓ H Statements H226
- ✓ P Statements P280

Sales unit :

40mL HYDRANAL™ - Water Standard 1.0
