


Specification

Phenol water-saturated, stabilized

A1624

Physical Description:	Liquid
Product Code:	A1624
Product Name:	Phenol water-saturated, stabilized
Short Description:	Stabilization: stabilized with 0.1 % 8-Hydroxyquinoline
Specifications:	Assay (titr.): approx. 75 % pH (20°C): 4.0 Heavy metals: max. 0.0005 % Water (K.F.): approx. 25 % Stability: approx. 15 months
Hazard pictograms	
UN:	2821
Class/PG:	6.1/II
ADR:	6.1/II
IMDG:	6.1/II
IATA:	6.1/II
WGK:	2
Storage:	2 - 8°C protected from light under argon

Signal Word:

Danger

Specification

Phenol water-saturated, stabilized

A1624

	GHS06 GHS08 GHS09
H Phrases:	EUH208 H302 H311+H331 H314 H341 H373 H411
P Phrases:	P280 P303+P361+P353 P305+P351+P338 P310 P321 P361+P364 P405 P501
CAS:	108-95-2
EINECS:	203-632-7
CS:	29071100
Index Nr.:	604-001-00-2

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340