

Specification

DNase I

A3778

Physical Description:	Solid
Product Code:	A3778
Product Name:	DNase I
Short Description:	delivery form: salt-free, freeze-dried powder, chromatographically prepared
Specifications:	Activity (Kunitz): min. 3000 U/mg
WGK:	1
Storage:	-20°C
Origin:	from bovine pancreas
M:	~31000 g/mol
CAS:	9003-98-9
EINECS:	232-667-0
CS:	35079090

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340