

Specification

EGTA for molecular biology

A0878

Physical Description:	Solid
Product Code:	A0878
Product Name:	EGTA for molecular biology
Specifications:	DNases/RNases/Proteases: not detectable Assay (titr.): min. 99 % Residue on ignition: max. 0.1 % Heavy metals (as Pb): max. 0.001 % Water: max. 1.0 % Fe: max. 0.001 %
WGK:	1
Storage:	RT
Molecular Formula:	$C_{14}H_{24}N_2O_{10}$
M:	380.35 g/mol
CAS:	67-42-5
EINECS:	200-651-2
CS:	29225000

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340