

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	34847-40ML-HON
<b>Description</b>	HYDRANAL™ - Water Standard 0.1

---

## Technical specifications :

- ✓ hydranal h2o stand. 0,1, stand. for k. fisch. 0,1mg/g, 4ml ampoule x10
- ✓ CAS N°95-47-6
- ✓ Molecular weight 106.17 g/mol
- ✓ Molecular formula C8H10
- ✓ Density 0.88 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point 1.013 hPa.32 °C
- ✓ Material form liquid
  
- ✓ Storage condition 2-8°C
- ✓ Package type Ampoule
- ✓ Package size 40mL
  
- ✓ GHS pictograms GHS07 ; GHS08 ; GHS02
- ✓ H Statements H412 ; H226 ; H304 ; H312 ; H332 ; H315 ; H319 ; H335
- ✓ P Statements P284 ; P280 ; P260 ; P210

## Sales unit :

40mL HYDRANAL™ - Water Standard 0.1

---