

Specification

Glycerol anhydrous for molecular biology

A2926

Refractive Index:	n ₂₀ /D 1.4740
Physical Description:	Liquid
Product Code:	A2926
Product Name:	Glycerol anhydrous for molecular biology
Specifications:	DNases/RNases/Proteases: not detectable Assay (GC): min. 99.5 % Acidity/Alkalinity: passes test Organic chlorine: max. 0.0005 % Fatty acids: max. 0.02 % Water (K.F.): max. 0.5 % Chloride: max. 0.0001 % Sulfate: max. 0.001 % As: max. 0.0001 % Fe: max. 0.0005 % Pb: max. 0.0001 %
WGK:	1
Storage:	RT
Origin:	from plants (non-animal origin!)
Molecular Formula:	C ₃ H ₈ O ₃
M:	92.10 g/mol
CAS:	56-81-5
EINECS:	200-289-5
CS:	29054500

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340