

TECHNICAL DATA SHEET

Part Number	34849-80ML-HON
Description	HYDRANAL™ - Water Standard 10.0

Technical specifications :

- ✓ hydranal h2o stand. 10.0, stand. for k. fisch. 10mg/g, 8ml ampoule x10
- ✓ CAS N°108-32-7
- ✓ Density 1.000 g/cm³ (20 °C)
- ✓ Flash point 32 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Ampoule
- ✓ Package size 80mL

- ✓ GHS pictograms GHS07 ; GHS05 ; GHS08 ; GHS02
- ✓ H Statements H318 ; H315 ; H304 ; H335 ; H336 ; H412 ; H226
- ✓ P Statements P260 ; P280 ; P284

Sales unit :

80mL HYDRANAL™ - Water Standard 10.0
