


Specification

Proteinase K solution 20 mg/ml

A4392

Physical Description:	Liquid
Product Code:	A4392
Product Name:	Proteinase K solution 20 mg/ml
Specifications:	DNases/RNases: not detectable Proteinase K: min. 600 mAnsonU/ml Composition: Glycerol (v/v): 40 % Proteinase K: 20 mg/mL CaCl ₂ · 2H ₂ O: 1 mM Tris · HCl (pH 7.5): 10 mM
Hazard pictograms	
WGK:	1
Storage:	2 - 8°C
Shipment:	wet ice
Signal Word:	Danger
GHS Symbols:	GHS08
H Phrases:	H334
P Phrases:	P261 P284 P304+P340 P342+P311 P501
CS:	38221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340