

Specification

TAE buffer (50X) for molecular biology

A4686

Physical Description:	Liquid
Product Code:	A4686
Product Name:	TAE buffer (50X) for molecular biology
Specifications:	DNases/RNases/Proteases: not detectable pH (20°C; H ₂ O): 8.5 ± 0.2 Composition: EDTA · Na ₂ · 2H ₂ O: 18.6 g/L (0,05 M) Acetic acid glacial: 60.05 g/L (1 M) Tris: 242.30 g/L (2 M)
WGK:	1
Storage:	RT
CS:	38221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340